

Product Description SALSA[®] MLPA[®] Probemix P295-B3 SPRED1

To be used with the MLPA General Protocol.

Version B3. As compared to version B2, two reference probes have been replaced. For complete product history see page 6.

Catalogue numbers:

- **P295-025R:** SALSA MLPA Probemix P295 SPRED1, 25 reactions.
- **P295-050R:** SALSA MLPA Probemix P295 SPRED1, 50 reactions.
- **P295-100R:** SALSA MLPA Probemix P295 SPRED1, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit, available for various number of reactions. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman capillary sequencers, respectively (see www.mlpa.com).

Certificate of Analysis: Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mlpa.com.

Precautions and warnings: For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mlpa.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

General Information: The SALSA MLPA Probemix P295 SPRED1 is a **research use only (RUO)** assay for the detection of deletions or duplications in the *SPRED1* gene, which is associated with Legius syndrome.

Legius syndrome (also known as Neurofibromatosis Type 1-like syndrome; NFLS) is characterised by multiple café-au-lait spots. Other features may include an abnormally large head (macrocephaly) and unusual facial characteristics. The syndrome has many symptomatic similarities with Neurofibromatosis Type I (NF1). *SPRED1* (Sprouty-related EVH1 domain containing 1) is a member of the Sprouty family of proteins that regulate growth factor-induced activation of the Ras/MAP kinase cascade. Mutations in the *SPRED1* gene cause Legius syndrome.

More information is available at <https://www.ncbi.nlm.nih.gov/books/NBK47312/>.

This SALSA MLPA Probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.

Gene structure and Transcript variants:

Entrez Gene shows transcript variants of each gene: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene>

For NM_ mRNA reference sequences: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide>

Locus Reference Genomic (LRG) database: <http://www.lrg-sequence.org/>

Probemix content: The SALSA MLPA Probemix P295-B3 SPRED1 contains 38 MLPA probes with amplification products between 149 and 454 nt. This includes 23 probes for the *SPRED1* gene. Multiple probes for each of the seven exons (or probes in close proximity of the exon) of the *SPRED1* gene are included. Probes in the upstream region and several probes in intronic regions of the gene have also been included in order to facilitate the detection of breakpoints of intragenic rearrangements. In addition, 15 reference probes are included and detect 15 different autosomal chromosomal locations. Complete probe sequences and the identity of the genes detected by the reference probes is available online (www.mlpa.com).

This Probemix contains nine quality control fragments generating amplification products between 64 and 121 nt: four DNA Quantity Fragments (Q-fragments), two DNA Denaturation Fragments (D-fragments), one benchmark fragment, one chromosome X, and one chromosome Y-specific fragment (see table below). More

information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mlpa.com.

Length (nt)	Name
64-70-76-82	Q-fragments (Only visible with <100 ng sample DNA)
88-96	D-fragments (Low signal of 88 or 96 fragment indicates incomplete denaturation)
92	Benchmark fragment
100	X-fragment (X chromosome specific)
105	Y-fragment (Y chromosome specific)

No DNA controls results in only five major peaks shorter than 121 nucleotides (nt): four Q-fragments at 64, 70, 76 and 82 nt, and one 19 nt peak corresponding to the unused portion of the fluorescent PCR primer. Non-specific peaks longer than 121 nt AND with a height >25% of the median of the four Q-fragments should not be observed. Note: peaks below this 25% threshold are not expected to affect MLPA reactions when sufficient amount of sample DNA (50-200 ng) is used.

MLPA technique: The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mlpa.com).

Required specimens: Extracted DNA, free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

Reference samples: All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from unrelated individuals who are from families without a history of Legius syndrome. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol.

Positive control DNA samples: MRC-Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Biobank (<https://catalog.coriell.org>) and DSMZ (<https://www.dsmz.de/home.html>) have a diverse collection of biological resources which may be used as a positive control DNA sample in your MLPA experiments. The quality of cell lines can change, therefore samples should be validated before use.

Data analysis: Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mlpa.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results: The standard deviation of all probes in the reference samples should be <0.10 and the dosage quotient (DQ) of the reference probes in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the DQ of the probes can be used to interpret MLPA results for autosomal or pseudo-autosomal chromosomes:

Copy Number status	Dosage quotient
Normal	$0.80 < DQ < 1.20$
Homozygous deletion	$DQ = 0$
Heterozygous deletion	$0.40 < DQ < 0.65$
Heterozygous duplication	$1.30 < DQ < 1.65$
Heterozygous triplication/ Homozygous duplication	$1.75 < DQ < 2.15$
Ambiguous copy number	All other values

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Incomplete DNA denaturation (e.g. due to salt contamination) can lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes are unlikely to have any relation to the condition tested for.

Limitations of the procedure:

- In most populations, the major cause of genetic defects in the *SPRED1* gene are small (point) mutations, most of which will not be detected by using SALSA MLPA Probemix P295 SPRED1.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNPs, point mutations, small indels) in the target sequence detected by a probe can cause false positive results. Mutations/SNPs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.

Confirmation of results: Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

SPRED1 mutation database: https://grenada.lumc.nl/LOVD2/mendelian_genes/home.php?select_db=SPRED1. We strongly encourage users to deposit positive results in the Leiden Open Variation Database (LOVD). Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report copy number changes detected by the reference probes, false positive results due to SNPs and unusual results (e.g., a duplication of *SPRED1* exons 5 and 7 but not exon 6) to MRC-Holland: info@mlpa.com.

Table 1. SALSA MLPA Probemix P295-B3 SPRED1

Length (nt)	SALSA MLPA probe	Chromosomal position (hg18) Reference <i>SPRED1</i>	
64-105	Control fragments – see table in probemix content section for more information		
149	Reference probe 10129-L11734	18q11	
154	SPRED1 probe 13044-L14227		Exon 3
160 Ø	SPRED1 probe 13045-L14228		Intron 4
166	Reference probe 08421-L08304	11q22	
172	SPRED1 probe 13046-L14229		Exon 5
178 Ø	SPRED1 probe 13047-L14230		Intron 2
185 Ø	SPRED1 probe 13048-L14231		Intron 1
196	SPRED1 probe 11017-L11686		Exon 4
202 *	Reference probe 03709-L25988	9q22	
208	SPRED1 probe 13049-L28437		Exon 4
214	SPRED1 probe 11015-L28438		Exon 2
220	Reference probe 12427-L13428	22q12	
230	SPRED1 probe 13050-L14233		Exon 1
238	SPRED1 probe 11019-L11688		Exon 6
247	SPRED1 probe 11016-L11685		Exon 3
257	Reference probe 17408-L21394	3p21	
265	SPRED1 probe 11013-L13179		Upstream
274 Ø	SPRED1 probe 13051-L14234		Intron 1
283	Reference probe 08680-L08692	13q32	
292	SPRED1 probe 13052-L14235		Exon 6
301	SPRED1 probe 11020-L11689		Exon 7
310 *	Reference probe 12442-L13443	14q24	
319	Reference probe 08048-L07829	5p15	
328	SPRED1 probe 13053-L14236		Exon 1
337	Reference probe 09937-L12248	8q13	
346	SPRED1 probe 11018-L11687		Exon 5
355	SPRED1 probe 13054-L14237		Exon 2
364	Reference probe 12989-L14146	19q13	
373	Reference probe 10718-L11300	6p12	
382	SPRED1 probe 13055-L14238		Exon 7
391	Reference probe 14043-L15641	7q36	
402	SPRED1 probe 11021-L11690		Exon 7
409	SPRED1 probe 13056-L14239		Upstream
420	Reference probe 08839-L13359	2p13	
427 Ø	SPRED1 probe 13057-L14240		Intron 5
436 Ø	SPRED1 probe 13058-L14241		Intron 2
445	Reference probe 12526-L13576	4q25	
454	Reference probe 08579-L08580	17q23	

* New in version B3 (from lot B3-0818 onwards).

Ø Intron probe. Only included to help determine the extent of a deletion/duplication. Copy number alterations of a single intron probe are unlikely to be related to the condition tested.

Note: The exon numbering used in this P295-B3 SPRED1 product description is the exon numbering from the RefSeq transcript NM_152594.2, which is identical to the NG_008980.1 sequence. The exon numbering and NM sequence used is from 09/2018, but can be changed (e.g. by NCBI) after the release of the product description. Please notify us of any mistakes: info@mlpa.com.

Table 2. SPRED1 probes arranged according to chromosomal location

Length (nt)	SALSA MLPA probe	SPRED1 exon	Ligation site NM_152594.2	Partial sequence (24 nt adjacent to ligation site)	Distance to next probe
		<i>start codon</i>	336-338 (Exon 1)		
409	13056-L14239	Upstream	839 nt before exon 1	GTGGCTCGACAA-CGATTGGGAAA	0.3 kb
265	11013-L13179	Upstream	509 nt before exon 1	CGGGCGTGAGGA-GCTGGCCCGCTG	0.5 kb
328	13053-L14236	Exon 1	49-50	GCCCCTCTCTTT-TTCCTTTCCAC	0.5 kb
230	13050-L14233	Exon 1	138 nt after exon 1	AGATAGCTGATT-TCTTTAGCCATT	8.4 kb
185 Ø	13048-L14231	Intron 1	8515 nt after exon 1	GGAATGTGGAAT-GTAATTCTTTTT	31.7 kb
274 Ø	13051-L14234	Intron 1	5928 nt before exon 2	TGGAACCTAGAT-CTGACTTGTAT	6.0 kb
214	11015-L28438	Exon 2	416-417	CGAGATGACTCA-AGTGGTGGATGG	0.1 kb
355	13054-L14237	Exon 2	509-510	TGTGCTGACTTT-TTTATCCGTGGA	6.7 kb
436 Ø	13058-L14241	Intron 2	6633 nt after exon 2	CTGTCTGCCTTT-TTCTCAAAGGT	12.2 kb
178 Ø	13047-L14230	Intron 2	3862 nt before exon 3	GTATTTATGTGA-GTCAAAGAGCTA	3.9 kb
154	13044-L14227	Exon 3	548-549	TGTCAGGTGGTT-TTGGAAATGTATG	0.1 kb
247	11016-L11685	Exon 3	663-664	CTGCTGATGCTA-GGGCTTTTGATA	2.3 kb
208	13049-L28437	Exon 4	117 nt before exon 4	TTGCCAGGCAGT-CCAGAAAGATCT	0.2 kb
196	11017-L11686	Exon 4	17 nt after exon 4	AATGGCTTGGA-GGAATTTGTAAA	10.6 kb
160 Ø	13045-L14228	Intron 4	4307 nt before exon 5	AAGACAGTGAAT-AGTAGGTATTGA	4.3 kb
172	13046-L14229	Exon 5	776-777	GAAGAGGATTCT-TCCAGTTCTCTA	0.1 kb
346	11018-L11687	Exon 5	839-840	AGTGAGCCTTAT-AGAAGCTCAAAT	6.8 kb
427 Ø	13057-L14240	Intron 5	2850 nt before exon 6	TTCATCCTAGAG-GTAACTATTACC	2.9 kb
238	11019-L11688	Exon 6	964-965	CAGAAGTATGGA-ATACGTACAGCG	0.1 kb
292	13052-L14235	Exon 6	84 nt after exon 6	TAATTCTCCATA-TAGTTGAATTGT	1.6 kb
301	11020-L11689	Exon 7	1174-1175	TGATTCCAGTAT-TCAGTTTCTAA	0.4 kb
402	11021-L11690	Exon 7	1570-1571	AGCTTTGTCTTT-CATTGTACCATG	0.2 kb
382	13055-L14238	Exon 7	1749-1750	AGCTTTTGGCAA-GCAATATGGAAT	
		<i>stop codon</i>	1668-1670 (Exon 7)		

Ø Intron probe. Only included to help determine the extent of a deletion/duplication. Copy number alterations of a single intron probe are unlikely to be related to the condition tested.

Note: The exon numbering used in this P295-B3 SPRED1 product description is the exon numbering from the RefSeq transcript NM_152594.2, which is identical to the NG_008980.1 sequence. The exon numbering and NM sequence used is from 09/2018, but can be changed (e.g. by NCB) after the release of the product description. Please notify us of any mistakes: info@mlpa.com.

Related SALSA MLPA probemixes

P044 NF2	Contains probes for the <i>NF2</i> gene, involved in Neurofibromatosis type 2.
P081/P082 NF1	Contains probes for the <i>NF1</i> gene, involved in Neurofibromatosis type 1.
P122 NF1 area	Contains probes for the 17q11.2 region, involved in Neurofibromatosis type 1.
P250 DiGeorge	Contains probes for the 22q11 region, used for primary screening of this region.
P258 SMARCB1	Contains probes for the <i>SMARCB1</i> gene. Deletions and mutations in <i>SMARCB1</i> are associated with malignant rhabdoid tumours and schwannomatosis.
P455 LZTR1	Contains probes for the <i>LZTR1</i> gene, involved in Neurofibromatosis type 3.

References

- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.

Selected publications using SALSA MLPA Probemix P295 SPRED1

- Evans DG et al. (2016). Comprehensive RNA analysis of the NF1 gene in classically affected NF1 affected individuals meeting NIH criteria has high sensitivity and mutation negative testing is reassuring in isolated cases with pigmentary features only. *EBioMedicine*. 7:212-220.
- Palma Milla C et al. (2018). Neurofibromatosis type I: mutation spectrum of NF1 in spanish patients. *Ann Hum Genet*.
- Pasmant E et al. (2015). Neurofibromatosis type 1 molecular diagnosis: what can NGS do for you when you have a large gene with loss of function mutations? *Eur J Hum Genet*. 23(5):596-601.
- Pinna V et al. (2015). p. Arg1809Cys substitution in neurofibromin is associated with a distinctive NF1 phenotype without neurofibromas. *Eur J Hum Genet*. 23(8):1068-1071.
- Van Minkelen R et al. (2014). A clinical and genetic overview of 18 years neurofibromatosis type 1 molecular diagnostics in the Netherlands. *Clin Genet*. 85(4):318-327.
- Zhang J et al. (2016). Identification of a PTPN11 hot spot mutation in a child with atypical LEOPARD syndrome. *Mol Med Rep*. 14(3):2639-2643.

P295 Product history	
Version	Modification
B3	Two reference probes have been replaced.
B2	One reference probe was added.
B1	One reference probe has been removed and six have been replaced, furthermore the control fragments have been adjusted (QDX2).
A1	First release.

Implemented changes in the product description

Version B3-01 – 03 October 2018 (01P)

- Product description adapted to a new product version (version number changed, changes in Table 1 and Table 2).
- Product description restructured and adapted to a new template.
- Small changes of probe lengths in Table 1 and 2 in order to better reflect the true lengths of the amplification products.

Version 09 – 02 December 2016 (55)

- Warning added below table 1 that an additional peak can appear at approximately 218 nt when the ligation reaction is at lower temperatures than recommended.

Version 08 – 30 April 2015 (54)

- Product description adapted to a new product version (version number changed, lot number added, small changes in Table 1 and Table 2, new pictures included).
- Minor textual changes.


Version 07 (49)

- Product description adapted to a new product version (version number changed, lot number added, small changes in Table 1 and Table 2, new pictures included).
- Minor textual changes.

Version 06 (48)

- Electropherogram pictures using the new MLPA buffer (introduced in December 2012) added.

More information: www.mlpa.com; www.mlpa.eu

	MRC-Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands
E-mail	info@mlpa.com (information & technical questions); order@mlpa.com (orders)
Phone	+31 888 657 200