

Product Description SALSA® MLPA® Probemix P003-D1 MLH1/MSH2

To be used with the MLPA General Protocol.

Version D1. As compared to version C1, two probes have been included that detect a recurrent 10 Mb inversion disrupting *MSH2*. In addition, one reference probe has been replaced (Table 1). For complete product history see page 12.

Catalogue numbers:

- **P003-025R:** SALSA MLPA probemix P003 MLH1/MSH2, 25 reactions.
- **P003-050R:** SALSA MLPA probemix P003 MLH1/MSH2, 50 reactions.
- **P003-100R:** SALSA MLPA probemix P003 MLH1/MSH2, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit, available for various number of reactions. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman capillary sequencers, respectively (see www.mlpa.com).

Certificate of Analysis: Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mlpa.com.

Precautions and warnings: For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mlpa.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

Intended use: This SALSA MLPA probemix P003 MLH1/MSH2 is an in vitro diagnostic (IVD)¹ or research use only (RUO) assay for the detection of exon deletion(s) or duplication(s) in specific regions of the *MLH1*, *MSH2* and *EPCAM* genes, as well as a recurrent 10 Mb inversion on chromosome arm 2p which disrupts the *MSH2* gene, in order to confirm a potential cause and clinical diagnosis for Lynch syndrome. This product can also be used for molecular genetic testing of at-risk family members.

This assay is for use with peripheral blood derived genomic DNA and not for use with DNA extracted from formalin-fixed paraffin embedded or fresh tumour materials. Deletions or duplications detected with the P003 MLH1/MSH2 probemix should be verified by using the SALSA MLPA probemix P248 MLH1/MSH2 Confirmation or a different technique. This verification must always be performed for copy number changes detected by only a single probe. P248 MLH1/MSH2 Confirmation cannot be used to verify deletions or duplications in *EPCAM*.

This assay is not intended to be used as a standalone assay for clinical decisions. Most defects in the *MLH1* and *MSH2* genes are point mutations, the majority of which will not be detected by MLPA. It is therefore recommended to use this SALSA MLPA probemix in combination with sequence analysis. The results of this test should be interpreted by a clinical molecular geneticist or equivalent.

¹Please note that this probemix is for In Vitro Diagnostic use (IVD) in the countries specified at the end of this product description. In all other countries, the product is for Research Use Only (RUO).

Clinical background: Germline defects in the *MLH1* and *MSH2* genes are the most frequent cause of a hereditary predisposition to Lynch syndrome (formerly known as hereditary non-polyposis colorectal cancer; HNPCC). Lynch syndrome is an inherited cancer of the digestive tract, particularly the colon and rectum. In addition, defects in the *MLH1* and *MSH2* genes increase the risk for numerous other cancers due to impaired DNA repair. More information is available on <http://www.ncbi.nlm.nih.gov/books/NBK1211/>.

Approximately 50% of Lynch syndrome cases are attributed to mutations in *MLH1*, 40% are attributed to mutations in *MSH2*, and 1-3% are due to *EPCAM* deletions. In addition, 10-20% of the cases can be explained by mutations in the *MSH6* and *PMS2* genes. Among the various defects in the *MLH1* and *MSH2* genes that have been found in patients are deletions and duplications of complete exons, which are usually

missed by standard sequence analysis. The MLPA technique can detect most of these deletions and duplications and therefore complements sequence analysis of the *MLH1* and *MSH2* genes.

Gene structure: The *MLH1* gene spans 57 kb on chromosome 3p22.2, and *MSH2* spans 80 kb on chromosome 2p21. The *MLH1* LRG_216 is identical to GenBank NG_007109.2 and the *MSH2* LRG_218 is identical to GenBank NG_007110.2. Both are available at www.lrg-sequence.org/.

Transcript variants:

MLH1 - Multiple transcript variants have been described: <http://www.ncbi.nlm.nih.gov/gene/4292>. *MLH1* transcript variant 1 (NM_000249.3, 2662 nt, coding sequence 199-2469) represents the most abundant transcript, encoding the full-length protein. This sequence is a reference standard in the NCBI RefSeqGene project. In Table 2a, the ligation sites of the *MLH1* MLPA probes are indicated according to this sequence which contains 19 exons. The ATG translation start site is located in exon 1 and the stop codon is located in exon 19.

MSH2 - Multiple transcript variants have been described: <http://www.ncbi.nlm.nih.gov/gene/4436>. *MSH2* transcript variant 1 (NM_000251.2, 3226 nt, coding sequence 126-2930) represents the longer transcript, encoding the longer isoform. This sequence is a reference standard in the NCBI RefSeqGene project. The ligation sites of the *MSH2* MLPA probes are indicated according to this sequence which contains 16 exons. The ATG translation start site is located in exon 1 and the stop codon is located in exon 16. The other transcript variant, NM_001258281.1, lacks an alternate segment of the first exon, including the translation start site, and is shorter at the N-terminus.

Exon numbering: The exon numbering used in this P003-D1 MLH1/MSH2 product description is the exon numbering for *MLH1* from the RefSeq transcript NM_000249.3, which is identical to the LRG_216 sequence, and for *MSH2* from the RefSeq transcript NM_000251.2, which is identical to the LRG_218 sequence. The exon numbering and NM sequence used is from 07/2018, but can be changed by NCBI after the release of the product description.

Probemix content: This SALSA MLPA probemix P003 MLH1/MSH2 version D1 contains 50 MLPA probes with amplification products between 130 and 499 nt (Table 1) including 39 probes for the *MLH1/MSH2/EPCAM* gene regions (Table 2), 2 probes that only generate signals on samples containing the recurrent 10 Mb inversion that disrupts *MSH2* (Wagner et al. 2002, Chen 2008, Rhees et al. 2014), and 9 reference probes that detect sequences outside this region. The identity of the genes detected by the reference probes is available online (www.mlpa.com).

One probe is present for each of the 19 exons of the *MLH1* gene and for each of the 16 exons of the *MSH2* gene. Two probes are included for exon 9 of *EPCAM* (formerly known as *TACSTD1*), a gene located just upstream of *MSH2*. Deletions of this last exon of *EPCAM* can result in silencing of the *MSH2* gene.

This Probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity Fragments (Q-fragments), two DNA Denaturation Fragments (D-fragments), one benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mlpa.com.

Length (nt)	Name
64-70-76-82	Q-fragments (Only visible with <100 ng sample DNA)
88-96	D-fragments (Low signal of 88 or 96 fragment indicates incomplete denaturation)
92	Benchmark fragment
100	X-fragment (X chromosome specific)
105	Y-fragment (Y chromosome specific)

MLPA technique: The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mlpa.com).

MLPA technique validation: Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation <0.10 for all probes over the experiment.

Required specimens: Extracted DNA from human peripheral blood, free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

Reference samples: All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from unrelated individuals who are from families without a history of hereditary predisposition to cancer. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol.

Positive control DNA samples: The NIBSC Institute in the U.K. provides a kit with 5 DNA samples containing *MLH1* or *MSH2* exon deletions or amplifications: <http://www.nibsc.org/documents/ifu/11-218-xxx.pdf>.

SALSA Binning DNA SD052: The SD052 Binning DNA provided with this probemix can be used as Binning DNA sample for binning of two inversion specific probes (265 nt probe 20091-SP0917-L28216 and 317 nt probe 20090-SP0916-L28222 for a 10 Mb inversion with breakpoint in intron 7 of *MSH2*). SD052 Binning DNA is a mixture of genomic DNA from healthy individuals and plasmid DNA that contains the target sequences detected by the above mentioned probes. Inclusion of one reaction with 5 μ l SD052 Binning DNA in initial MLPA experiments is essential as it can be used to aid in data binning of the peak pattern using Coffalyser.Net software. Furthermore, Binning DNA should be included in the experiment whenever changes have been applied to the set-up of the capillary electrophoresis device (e.g. when capillaries have been renewed). Binning DNA should never be used as a reference sample in the MLPA data analysis, neither should it be used in quantification of mutation signals, as for this purpose true inversion positive patient samples or cell lines should be used. It is strongly advised to use DNA sample and reference DNA samples extracted with the same method and derived from the same source of tissue. For further details, please consult the SD052 Binning DNA product description provided.

Performance characteristics: Large deletions account for ~20% of variants in *MSH2*, ~5% of variants in *MLH1* and for all *EPCAM* variants that cause Lynch syndrome (<http://www.ncbi.nlm.nih.gov/books/NBK1211/>), and the great majority of these can be detected by this P003 MLPA probemix. Furthermore, the recurrent 10 Mb inversion of chromosome 2 which has one breakpoint in intron 7 of the *MSH2* gene can also be detected by this probemix (Wagner et al. 2002, Chen 2008, Rhee et al. 2014). *MLH1* or *MSH2* duplications have been rarely reported, but should also be detectable with this probemix. Note that the P003 *MLH1/MSH2* probemix does not contain probes for the detection of *MSH6* and *PMS2* genes, which are also implicated in Lynch syndrome. The analytical sensitivity and specificity for the detection of deletions in the aforementioned genes using point mutation negative samples with Lynch syndrome is very high and can be considered $>99\%$.

Analytical performance can be compromised by: SNPs or other polymorphisms (e.g. indels) in the DNA target sequence, impurities in the DNA sample, incomplete DNA denaturation, the use of insufficient or too much sample DNA, the use of insufficient or unsuitable reference samples, problems with capillary electrophoresis or a poor data normalisation procedure and other technical errors. The MLPA General Protocol contains technical guidelines and information on data evaluation/normalisation.

Data analysis: Coffalyser.Net software must be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mlpa.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results: The expected results for MLH1/MSH2 region specific MLPA probes are allele copy numbers of 2 (normal), 1 (heterozygous deletion), 3 (heterozygous duplication) and occasionally 4 (heterozygous triplication). Homozygous deletions may occur, but are extremely rare.

The standard deviation of all probes in the reference samples should be <0.10 and the dosage quotient (DQ) of the reference probes in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the DQ of the probes can be used to interpret MLPA results for autosomal or pseudo-autosomal chromosomes:

Copy Number status	Dosage quotient
Normal	$0.80 < DQ < 1.20$
Homozygous deletion	$DQ = 0$
Heterozygous deletion	$0.40 < DQ < 0.65$
Heterozygous duplication	$1.30 < DQ < 1.65$
Heterozygous triplication/ Homozygous duplication	$1.75 < DQ < 2.15$
Ambiguous copy number	All other values

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Incomplete DNA denaturation (e.g. due to salt contamination) can lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might in some cases not result in inactivation of that gene copy.
- Copy number changes detected by reference probes are unlikely to have any relation to the condition tested for.

Notes P003 MLH1-MSH2 results:

- Lynch syndrome due to *MLH1/MSH2/EPCAM* gene defects is an autosomal dominant disorder. Inactivation of a single copy of the *MLH1/MSH2* gene or deletion of the last *EPCAM* exon is expected to be pathogenic.
- A heterozygous deletion of one or more *MLH1* or *MSH2* exons that are present in the major transcript variants, NM_000249.3 (*MLH1*) and NM_000251.2 (*MSH2*), is expected to result in Lynch syndrome.
- Simultaneous deletion of the 472 and 481 nt *EPCAM* probes is a strong indication that the last *EPCAM* exon (exon 9) is disrupted, which can lead to methylation and inactivation of *MSH2* (Ligtenberg et al. 2009, Kovacs et al. 2009).
- The presence of peaks at 265 and 317 nt, or only at 265 nt, is an indication that the 10 Mb chromosome 2 inversion is present in the sample DNA (Wagner et al. 2002, Chen 2008, Rhees et al. 2014). These probes do not generate a signal on normal samples. Please note that the peak height of these two probes can be strongly affected by DNA fragmentation. The 20091-SP0917-L28216 probe present at 265 nt and the 20090-SP0916-L28222 probe present at 317 nt require an intact DNA fragment to bind to of at least 220 nt or 520 nt, respectively, in inversion positive DNA samples, which is much longer than the length requirements for other MLPA probes (60-80 nt). DNA samples in which both probes show a clear peak, as well as samples in which only a clear peak for one of the two probes is observed, should be further investigated even when these peaks are much lower as obtained on the SD052 Binning DNA

sample. A non-specific peak close to the 265 nt peak can appear if the ligation reaction is performed at a temperature below 54°C.

- A duplication of an internal part of a gene usually results in a defective copy of that gene, as the duplicated sequence is typically located directly adjacent to the original sequence, resulting in a defective transcript. However, duplication of the *complete MLH1* or *MSH2* gene is not expected to result in disease. Please note the remark above on duplications that include the first or last exon of a gene.

Limitations of the procedure:

- In most populations, the major cause of genetic defects in the *MLH1* and *MSH2* gene are small (point) mutations, most of which will not be detected by using SALSA MLPA probemix P003 MLH1/MSH2.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and does not detect translocations or copy number neutral inversions, except the 10 Mb inversion with breakpoint in intron 7 of *MSH2*. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in the genes or chromosomal regions *do* exist but remain undetected.
- Sequence changes (e.g. SNPs, point mutations, small indels) in the target sequence detected by a probe can cause false positive results. Mutations/SNPs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.

Confirmation of results: Deletions or duplications obtained with the P003 MLH1/MSH2 probemix must be verified by using SALSA MLPA probemix P248 MLH1/MSH2 Confirmation or a different technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH. Almost all probes included in probemix P248 MLH1/MSH2 Confirmation have different ligation sites from those in P003, except for the 148 nt exon 1 *MSH2* probe. The P072 MSH6 MLPA probemix contains additional probes for *EPCAM* and the region between *MSH2* and *EPCAM*. In addition, the P072 probemix can be used to detect copy number variations in the *MSH6* gene.

Copy number changes detected by only a single probe always require confirmation by the probemix P248 MLH1/MSH2 Confirmation or by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

MLH1/MSH2 mutation database: <http://www.insight-group.org/>; <http://www.lovd.nl>. We encourage users to deposit positive results in the LOVD database of the International Society for Gastrointestinal Hereditary Tumours (INSiGHT). Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report copy number changes detected by the reference probes, false positive results due to SNPs and unusual results (e.g., a duplication of *MLH1* exons 6 and 8 but not exon 7) to MRC-Holland: info@mlpa.com.

Table 1. SALSA MLPA Probemix P003-D1 MLH1/MSH2

Length (nt)	SALSA MLPA probe	Chromosomal position (hg18) ^(a)		
		reference	MLH1	MSH2 / EPCAM
64-105	Control fragments – see table in probemix content section for more information			
130	Reference probe 00797-L21056	5q31		
136	Reference probe 00981-L00566	10p11		
142	MLH1 probe 14701-L19726		Exon 1	
148 »	MSH2 probe 12036-L02162			Exon 1
154	MLH1 probe 01008-L00577		Exon 2	
160	MSH2 probe 00906-L00494			Exon 2
166	MLH1 probe 00888-L00476		Exon 3	
172	MSH2 probe 01029-L00601			Exon 3
178	MLH1 probe 00889-L23928		Exon 4	
184	MSH2 probe 00908-L00496			Exon 4
190	MSH2 probe 11287-L12006			Exon 1
196	Reference probe 17056-L20134	7q31		
203 ¥	MLH1 probe 20360-L28308		Exon 5	
211	MSH2 probe 00909-L00497			Exon 5
217	MSH2 probe 13145-L14624			Exon 1
224 ¥	MLH1 probe 14702-L28215		Exon 6	
231	MSH2 probe 15298-L17064			Exon 6
236	Reference probe 15941-L18067	2q24		
242	MLH1 probe 00892-L16355		Exon 7	
249 «	MSH2 probe 11634-L16356			Exon 7
256	MLH1 probe 00893-L00481		Exon 8	
265 * §	Inversion probe 20091-SP0917-L28216			Intron 7
269 ¥	MSH2 probe 00912-L28217			Exon 8
278 ¥ +	MLH1 probe 00894-L28218		Exon 9	
287 *	Reference probe 18920-L25191	1p36		
292 ¥	MSH2 probe 00913-L28221			Exon 9
301	MLH1 probe 00895-L00483		Exon 10	
310	MSH2 probe 11288-L12007			Exon 10
317 * §	Inversion probe 20090-SP0916-L28222			Intron 7
326 ¥	MLH1 probe 00896-L18364		Exon 11	
332 ¥	MSH2 probe 00915-L28223			Exon 11
340 ¥ ∫	MLH1 probe 14703-L28224		Exon 12	
346	MSH2 probe 18133-L23925			Exon 12
355	MLH1 probe 00898-L23926		Exon 13	
364	MSH2 probe 01013-L00575			Exon 13
375	Reference probe 00681-L11147	4q25		
382	MLH1 probe 00899-L00586		Exon 14	
391	MSH2 probe 00918-L00506			Exon 14
401	MLH1 probe 00900-L00488		Exon 15	
409	MSH2 probe 00919-L00585			Exon 15
418	MLH1 probe 01009-L00576		Exon 16	
427	MSH2 probe 01053-L14623			Exon 16
436	MLH1 probe 01030-L00602		Exon 17	
445	MLH1 probe 01031-L00603		Exon 18	
454	MLH1 probe 12094-L12994		Exon 19	
463	Reference probe 00979-L00568	10p14		
472	EPCAM probe 13147-L14404			(Exon 9)
481	EPCAM probe 18132-L24050			(Exon 9)
490 «	Reference probe 04274-L24051	13q12		
499	Reference probe 14882-L21050	14q11		

(a) The exon numbering used in this P003-D1 MLH1/MSH2 product description is the exon numbering for *MLH1* from the RefSeq transcript NM_000249.3, which is identical to the LRG_216 sequence, and for *MSH2* from the RefSeq transcript NM_000251.2, which is identical to the LRG_218 sequence. The exon numbering

and NM sequence used is from 07/2018, but can be changed by NCBI after the release of the product description.

* New in version D1 (from lot D1-0815 onwards).

‡ Changed in version D1 (from lot D1-0815 onwards). Small change in length, no change in sequence detected.

« Probe located within, or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

+ See note on the 278 nt *MLH1* exon 9 probe below Table 2a.

‡ See note on the 340 nt *MLH1* exon 12 probe below Table 2a.

§ Mutation-specific probe. This probe will only generate a signal when the 10 Mb inversion with breakpoint in intron 7 of *MSH2* is present. This probe consists of three parts and has two ligation sites. See note about this probe on page 4, section interpretation of results.

» Detects the same sequence as the 355 nt probe (02735-L20467) in SALSA MLPA probemix P248 *MLH1/MSH2* Confirmation.

Table 2. P003 probes arranged according to chromosomal location

Table 2a. *MLH1* gene

Length (nt)	SALSA MLPA probe	<i>MLH1</i> Exon ^(a)	Ligation site ^(b) NM_000249.3	Partial sequence ^(c) (24 nt adjacent to ligation site)	Distance to next probe
		<i>start codon</i>	<i>199-201 (exon 1)</i>		
142	14701-L19726	Exon 1	205-206	CCAAAATGTCGT-TCGTGGCAGGGG	3.1 kb
154	01008-L00577	Exon 2	383-384	GATTCAGATCCA-AGACAATGGCAC	4.3 kb
166	00888-L00476	Exon 3	463-464	TGCAGTCCTTTG-AGGATTTAGCCA	3.4 kb
178	00889-L23928	Exon 4	527-528	CATAAGCCATGT-GGCTCATGTTAC	2.6 kb
203	20360-L28308	Exon 5	627-628	AAACCATGTGCT-GGCAATCAAGGG	1.8 kb
224	14702-L28215	Exon 6	664-665	TGGAGGACCTTT-TTTACAACATAG	3.0 kb
242	00892-L16355	Exon 7	764-765	ACACAATGCAGG-CATTAGTTTCTC	0.2 kb
256	00893-L00481	Exon 8	818-819	TGTTAGGACACT-ACCCAATGCCTC	2.4 kb
278 +	00894-L28218	Exon 9	920-921	CCTAGCCTTCAA-AATGAATGGTTA	3.1 kb
301	00895-L00483	Exon 10	1019-1020	TTCTTGAGAAA-AGCCATAGAAAC	2.9 kb
326	00896-L18364	Exon 11	1190-1191	GCAGCACATCGA-GAGCAAGCTCCT	5.5 kb
340 ‡	14703-L28224	Exon 12	1482-1483	AGGCAGCAAGAT-GAGGAGATGCTT	3.0 kb
355	00898-L23926	Exon 13	1665-1666	TCCCGAAAGGAA-ATGACTGCAGCT	11.4 kb
382	00899-L00586	Exon 14	1799-1800	CGTGGGCTGTGT-GAATCCTCAGTG	2.1 kb
401	00900-L00488	Exon 15	1892-1893	CCAGATACTCAT-TTATGATTTTGC	5.3 kb
418	01009-L00576	Exon 16	1967-1968	CATGCTTGCTT-AGATAGTCCAGA	1.0 kb
436	01030-L00602	Exon 17	2172-2173	ATCTTCATTCTT-CGACTAGCCACT	0.4 kb
445	01031-L00603	Exon 18	2269-2268 reverse	CTCCTCAGATAT-GTACTGCTTCCG	1.8 kb
454	12094-L12994	Exon 19	2595-2594 reverse	TATCAGAAGGCA-AGTATAAGTCTT	
		<i>stop codon</i>	<i>2467-2469 (exon 19)</i>		

(a) The exon numbering used in this P003-D1 *MLH1/MSH2* product description is the exon numbering for *MLH1* from the RefSeq transcript NM_000249.3, which is identical to the LRG_216 sequence. The exon numbering and NM sequence used is from 07/2018, but can be changed by NCBI after the release of the product description.

(b) Ligation sites of the P003 *MLH1* MLPA probes are indicated according to RefSeq sequence NM_000249.3 containing 19 exons.

(c) Only partial probe sequences are shown. Complete probe sequences are available at www.mlpa.com. Please notify us of any mistakes: info@mlpa.com.

+ The *MLH1* exon 9 probe in the previous version could give a duplication signal that could not be reproduced by the *MLH1* exon 9 probe in the P248 *MLH1/MSH2* Confirmation probemix. This apparent duplication was probably caused by the c.702G>A variant (rs35908749), located 20 nt from the ligation site. In the current D1 version (lot 0815 onwards) this probe has been adjusted, resulting in normal results when a heterozygous positive sample for this SNP is analysed. Please notify us if you still obtain different results with this probe than with the corresponding probe in the P248 probemix.

¶ The 340 nt MLH1 exon 12 probe has been found to give false positive deletion/duplication results. This is probably due to an unusual low sensitivity of this probe to sample DNA depurination. Depurination of sample DNA can occur in FFPE derived samples or in samples with insufficient buffering capacity. The MLH1 exon 12 probe in the P248 probemix does not have this problem and can be used for confirmation. More information on sample DNA depurination is available on www.mlpa.com.

Table 2b. *MSH2* gene

Length (nt)	SALSA MLPA probe	<i>MSH2</i> Exon ^(a)	Ligation ^(b) site	Partial sequence ^(c) (24 nt adjacent to ligation site)	Distance to next probe
<i>EPCAM</i> gene					
			NM_002354.2		
			<i>EPCAM</i> stop codon	<i>1301-1303 (exon 9)</i>	
481 #	18132-L24050	<i>EPCAM</i> Exon 9	1352-1353	AAATGGACACAA-ATTACAAATGTG	0.1 kb
472	13147-L14404	<i>EPCAM</i> Exon 9	1483-1482 reverse	GGTCAAATTTCA-AGATTGGTAAAG	16.0 kb to <i>MSH2</i>
<i>MSH2</i> gene					
			NM_000251.2		
		<i>start codon</i>	<i>126-128 (exon 1)</i>		
217	13145-L14624	Exon 1	257 nt before exon 1	CCGGGCACATTA-CGAGCTCAGTGC	0.2 kb
148 »	12036-L02162	Exon 1	67 nt before exon 1	GCGTGCGCGGGA-AGCTGGGCCGCG	0.7 kb
190	11287-L12006	Exon 1	269 nt after exon 1 reverse	GAAGTAGAACAA-TGCATTAATAATG	4.8 kb
160	00906-L00494	Exon 2	428-429	TATAGAGTTGAA-GTTTATAAGAAT	1.7 kb
172	01029-L00601	Exon 3	576-577	TTGTGGGTGTTA-AAATGTCCGCAG	2.3 kb
184	00908-L00496	Exon 4	863-864	CGGTTGTTGAAA-GGCAAAAAGGGA	1.8 kb
211	00909-L00497	Exon 5	943-944	ACTGTCTGCGGT-AATCAAGTTTTT	2.1 kb
231	15298-L17064	Exon 6	1118-1119	GCCTTGCTGAAT-AAGTGTAACACC	13.5 kb
249 «	11634-L16356	Exon 7	1322-1323	AGACAAGCAGCA-AACTTACAAGAT	15.7 kb
269	00912-L28217	Exon 8	1449-1450	CTCCTCTTACTG-ATCTTCGTTCTG	17.5 kb
292	00913-L28221	Exon 9	1536-1537	AATTCCTTGTA-AACCTTCATTTG	3.9 kb
310	11288-L12007	Exon 10	189 nt after exon 10	GACTGAAGTGGT-ACTTTGGGTCTA	4.0 kb
332	00915-L28223	Exon 11	1832-1831 reverse	GCTTCTTCATAT-TCTGTTTTATTT	4.1 kb
346	18133-L23925	Exon 12	1980-1981	CACCTGTTCCAT-ATGTACGACCAG	1.3 kb
364	01013-L00575	Exon 13	2196-2197	TCATGGCCCAA-TTGGGTGTTTTG	1.9 kb
391	00918-L00506	Exon 14	2403-2404	CCTACGATGGAT-TTGGGTAGCAT	2.5 kb
409	00919-L00585	Exon 15	2686-2687	ACTTGAGGAGTT-TCAGTATATTGG	2.0 kb
427	01053-L14623	Exon 16	2762-2763	GTGTTTCAGCAA-GGTGAAAAAATT	
		<i>stop codon</i>	<i>2928-2930 (exon 16)</i>		
265 §	20091-0917-L28216	Intron 7	This probe will only generate a signal on mutation positive DNA samples (10 Mb inversion with breakpoint in <i>MSH2</i> intron 7)		
317 §	20090-0916-L28222	Intron 7	This probe will only generate a signal on mutation positive DNA samples (10 Mb inversion with breakpoint in <i>MSH2</i> intron 7)		

(a) The exon numbering used in this P003-D1 MLH1/*MSH2* product description is the exon numbering for *MSH2* from the RefSeq transcript NM_000251.2, which is identical to the LRG_218 sequence. The exon numbering and NM sequence used is from 07/2018, but can be changed by NCBI after the release of the product description.

(b) Ligation sites of the P003 *MSH2* MLPA probes are indicated according to RefSeq sequence NM_000251.2 containing 16 exons.

(c) Only partial probe sequences are shown. Complete probe sequences are available at www.mlpa.com. Please notify us of any mistakes: info@mlpa.com.

« Probe located within, or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

§ Mutation-specific probe. This probe will only generate a signal when the 10 Mb inversion with breakpoint in intron 7 of *MSH2* is present. This probe consists of three parts and has two ligation sites. See note about this probe on page 4, section interpretation of results.

» Detects the same sequence as the 355 nt probe (02735-L20467) in SALSA MLPA probemix P248 MLH1/MSH2 Confirmation.

This probe's specificity relies on a single nucleotide difference compared to a related gene or pseudogene. As a result, an apparent duplication of only this probe can be the result of a non-significant single nucleotide sequence change in the related gene or pseudogene.

Related SALSA MLPA probemixes

COLORECTAL CANCER (CRC) PROBEMIXES	Condition		Gene	Probemix
	Lynch Syndrome (HNPCC)		<i>MLH1</i>	P003, ME011, P248 (Confirmation of P003), ME042
			<i>MSH2</i>	P003, ME011, P248 (Confirmation of P003)
			<i>MSH6</i>	P072, ME011
			<i>PMS2</i>	P008, ME011
			<i>EPCAM</i>	P003, P072
	Polyposis Syndrome	MAP	<i>MUTYH</i>	P378, P043, P072
		AFAP	<i>APC</i>	P043
		FAP	<i>APC</i>	P043

References

- Chen JM (2008). The 10-Mb paracentric inversion of chromosome arm 2p in activating *MSH2* and causing hereditary nonpolyposis colorectal cancer: re-annotation and mutational mechanisms. *Genes Chromosomes Cancer*. 47:543-5.
- Kovacs ME et al. (2009). Deletions removing the last exon of *TACSTD1* constitute a distinct class of mutations predisposing to Lynch syndrome. *Hum Mutat*. 30:197-203.
- Ligtenberg MJ et al. (2009). Heritable somatic methylation and inactivation of *MSH2* in families with Lynch syndrome due to deletion of the 3' exons of *TACSTD1*. *Nat Genet*. 41:112-7.
- Rhees J et al. (2014). Inversion of exons 1-7 of the *MSH2* gene is a frequent cause of unexplained Lynch syndrome in one local population. *Fam Cancer*. 13:219-25.
- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res*. 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat*. 28:205.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem*. 421:799-801.
- Wagner A et al. (2002). A 10-Mb paracentric inversion of chromosome arm 2p inactivates *MSH2* and is responsible for hereditary nonpolyposis colorectal cancer in a North-American kindred. *Genes Chromosomes Cancer*. 35:49-57.

Selected publications using SALSA® MLPA® Probemix P003 MLH1/MSH2

- Abu Freha N et al. (2018). Constitutional mismatch repair deficiency and Lynch syndrome among consecutive Arab Bedouins with colorectal cancer in Israel. *Fam Cancer*. 17:79-86.

- Akbari MR et al. (2017). Correlation between germline mutations in MMR genes and microsatellite instability in ovarian cancer specimens. *Fam Cancer*. 16:351-5.
- Alonso-Espinaco V et al. (2011). Novel MLH1 duplication identified in Colombian families with Lynch syndrome. *Genet Med*. 13:155-60.
- Baas AF et al. (2013). Agenesis of the corpus callosum and gray matter heterotopia in three patients with constitutional mismatch repair deficiency syndrome. *Eur J Hum Genet*. 21:55-61.
- Bashyam MD et al. (2015). Evidence for presence of mismatch repair gene expression positive Lynch syndrome cases in India. *Mol Carcinog*. 54:1807-14.
- Berzina D et al. (2012). Novel germline MLH1 and MSH2 mutations in Latvian Lynch syndrome families. *Exp Oncol*. 34:49-52.
- Borelli I et al. (2013). A unique MSH2 exon 8 deletion accounts for a major portion of all mismatch repair gene mutations in Lynch syndrome families of Sardinian origin. *Eur J Hum Genet*. 21:154-61.
- Borelli I et al. (2014). A founder MLH1 mutation in Lynch syndrome families from Piedmont, Italy, is associated with an increased risk of pancreatic tumours and diverse immunohistochemical patterns. *Fam Cancer*. 13:401-13.
- Carneiro da Silva F et al. (2015). Clinical and Molecular Characterization of Brazilian Patients Suspected to Have Lynch Syndrome. *PLoS One*. 10:e0139753.
- Chika N et al. (2016). Prevalence of Lynch syndrome and Lynch-like syndrome among patients with colorectal cancer in a Japanese hospital-based population. *Jpn J Clin Oncol*. 47:108-17.
- Cini G et al. (2015). Concomitant mutation and epimutation of the MLH1 gene in a Lynch syndrome family. *Carcinogenesis*. 36:452-8.
- Czink E et al. (2017). Successful immune checkpoint blockade in a patient with advanced stage microsatellite unstable biliary tract cancer. *Cold Spring Harb Mol Case Stud*. mcs. a001974.
- De Lellis L et al. (2013). Integrative analysis of hereditary nonpolyposis colorectal cancer: the contribution of allele-specific expression and other assays to diagnostic algorithms. *PLoS One*. 8:e81194.
- Dominguez-Valentin M et al. (2013). Mutation spectrum in South American Lynch syndrome families. *Hered Cancer Clin Pract*. 11:18.
- Duraturo F et al. (2013). Contribution of large genomic rearrangements in Italian Lynch syndrome patients: characterization of a novel alu-mediated deletion. *Biomed Res Int*. 2013:219897.
- Egoavil C et al. (2013). Prevalence of Lynch syndrome among patients with newly diagnosed endometrial cancers. *PLoS One*. 8:e79737.
- Eguchi H et al. (2016). Identification of a Japanese Lynch syndrome patient with large deletion in the 3' region of the EPCAM gene. *Jpn J Clin Oncol*. 46:178-84.
- Goldberg Y et al. (2014). Lynch Syndrome in high risk Ashkenazi Jews in Israel. *Fam Cancer*. 13:65-73.
- Huth C et al. (2012). The molecular basis of EPCAM expression loss in Lynch syndrome-associated tumors. *Mod Pathol*. 25:911-6.
- Kamiza AB et al. (2015). Risk factors associated with colorectal cancer in a subset of patients with mutations in MLH1 and MSH2 in Taiwan fulfilling the Amsterdam II criteria for Lynch syndrome. *PLoS one*. 10:e0130018.
- Kato A et al. (2016). Isolated Loss of PMS2 Immunohistochemical Expression is Frequently Caused by Heterogenous MLH1 Promoter Hypermethylation in Lynch Syndrome Screening for Endometrial Cancer Patients. *Am J Surg Pathol*. 40:770-6.
- Kidambi TD et al. (2016). Constitutional MLH1 methylation presenting with colonic polyposis syndrome and not Lynch syndrome. *Fam Cancer*. 15:275-80.
- Kloor M et al. (2011). Analysis of EPCAM protein expression in diagnostics of Lynch syndrome. *J Clin Oncol*. 29:223-7.
- Kovac MB et al. (2015). High-Resolution Breakpoint Analysis Provides Evidence for the Sequence-Directed Nature of Genome Rearrangements in Hereditary Disorders. *Human mutation*. 36:250-9.
- Kraus C et al. (2015). Comprehensive screening for mutations associated with colorectal cancer in unselected cases reveals penetrant and nonpenetrant mutations. *Int. J Cancer*. 136.
- Leclerc J et al. (2018). Diversity of genetic events associated with MLH1 promoter methylation in Lynch syndrome families with heritable constitutional epimutation. *Genet Med*.
- Liu Y et al. (2014). Systematic study on genetic and epimutational profile of a cohort of Amsterdam criteria-defined Lynch Syndrome in Singapore. *PLoS One*. 9:e94170.
- Loconte DC et al. (2014). A rare MSH2 mutation causes defective binding to hMSH6, normal hMSH2 staining, and loss of hMSH6 at advanced cancer stage. *Hum Pathol*. 45:2162-7.
- Loizidou MA et al. (2014). The Mutational Spectrum of Lynch Syndrome in Cyprus. *PLoS one*. 9:e105501.

- Maccaroni E et al. (2015). Prognostic impact of mismatch repair genes germline defects in colorectal cancer patients: are all mutations equal? *Oncotarget*. 6:38737-48.
- Magnani G et al. (2015). Molecular features and methylation status in early onset (≤ 40 years) colorectal cancer: a population based, case-control study. *Gastroenterol Res Pract*. 2015: 132190.
- Malesci A et al. (2014). Molecular heterogeneity and prognostic implications of synchronous advanced colorectal neoplasia. *Br J Cancer*. 110:1228-35.
- Maresca L et al. (2015). MSH2 role in BRCA1-driven tumorigenesis: A preliminary study in yeast and in human tumors from BRCA1-VUS carriers. *Eur J Med Genet*. 58:531-9.
- Mensenkamp AR et al. (2014). Somatic mutations in MLH1 and MSH2 are a frequent cause of mismatch-repair deficiency in Lynch syndrome-like tumors. *Gastroenterology*. 146:643-6 e648.
- Moir-Meyer GL et al. (2015). Rare germline copy number deletions of likely functional importance are implicated in endometrial cancer predisposition. *Human genetics*. 134:269.
- Mur P et al. (2014). Identification of a founder EPCAM deletion in Spanish Lynch syndrome families. *Clin Genet*. 85:260-6.
- Musulen E et al. (2013). Usefulness of epithelial cell adhesion molecule expression in the algorithmic approach to Lynch syndrome identification. *Hum Pathol*. 44:412-6.
- Perez-Cabornero L et al. (2011a). Characterization of new founder Alu-mediated rearrangements in MSH2 gene associated with a Lynch syndrome phenotype. *Cancer Prev Res (Phila)*. 4:1546-55.
- Perez-Cabornero L et al. (2011b). Frequency of rearrangements in Lynch syndrome cases associated with MSH2: characterization of a new deletion involving both EPCAM and the 5' part of MSH2. *Cancer Prev Res (Phila)*. 4:1556-62.
- Petersen SM et al. (2013). Functional examination of MLH1, MSH2, and MSH6 intronic mutations identified in Danish colorectal cancer patients. *BMC Med Genet*. 14:103.
- Pinheiro M et al. (2011). A novel exonic rearrangement affecting MLH1 and the contiguous LRRFIP2 is a founder mutation in Portuguese Lynch syndrome families. *Genet Med*. 13:895-902.
- Pinto C et al. (2016). Co-occurrence of nonsense mutations in MSH6 and MSH2 in Lynch syndrome families evidencing that not all truncating mutations are equal. *J Hum Genet*. 61:151-6.
- Rey J-M et al. (2017). Improving Mutation Screening in Patients with Colorectal Cancer Predisposition Using Next-Generation Sequencing. *J Mol Diagn*. 19:589-601.
- Rohlin A et al. (2016). GREM1 and POLE variants in hereditary colorectal cancer syndromes. *Genes, Chromosomes and Cancer*. 55:95-106.
- Romero A et al. (2013). Frequency and variability of genomic rearrangements on MSH2 in Spanish Lynch Syndrome families. *PLoS One*. 8:e72195.
- Rubio I et al. (2016). Analysis of Lynch Syndrome Mismatch Repair Genes in Women with Endometrial Cancer. *Oncology*. 91:171-6.
- Schneider NB et al. (2018). Germline MLH1, MSH2 and MSH6 variants in Brazilian patients with colorectal cancer and clinical features suggestive of Lynch Syndrome. *Cancer Med*. 7:2078-88.
- Shiozawa M et al. (2013). Partial duplication of MSH2 spanning exons 7 through 14 in Lynch syndrome. *J Gastroenterol*. 48:770-6.
- Sloane MA et al. (2015). Mosaic Epigenetic Inheritance as a Cause of Early-Onset Colorectal Cancer. *JAMA Oncol*. 1:953-7.
- Smith MJ et al. (2016). The contribution of whole gene deletions and large rearrangements to the mutation spectrum in inherited tumor predisposing syndromes. *Human mutation*. 37:250-256.
- Soares BL et al. (2018). Screening for germline mutations in mismatch repair genes in patients with Lynch syndrome by next generation sequencing. *Fam Cancer*. 17:387-94.
- Stigliano V et al. (2014). Early-onset colorectal cancer patients without family history are "at very low risk" for lynch syndrome. *J Exp Clin Cancer Res*. 33:1.
- Sun KK et al. (2016). Small bowel adenocarcinoma in Lynch syndrome: A case report. *Oncol Lett*. 12:1602-4.
- Suzuki O et al. (2017). Prevalence and clinicopathologic/molecular characteristics of mismatch repair-deficient colorectal cancer in the under-50-year-old Japanese population. *Surg Today*. 47:1135-46.
- Takahashi K et al. (2017). Clinical characteristics of Lynch-like cases collaterally classified by Lynch syndrome identification strategy using universal screening in endometrial cancer. *Gynecol Oncol*. 147:388-95.
- Tanyi M et al. (2014). MLH1 and MSH2 mutation screening in HNPCC families of Hungary - Two new MMR gene mutations. *Eur J Surg Oncol*. 40:1445-52.

- Urakami S et al. (2018). Clinicopathological characteristics of patients with upper urinary tract urothelial cancer with loss of immunohistochemical expression of the DNA mismatch repair proteins in universal screening. *Int J Urol*. 25:151-6.
- Wilders EA et al. (2014). Functional analysis of MSH2 unclassified variants found in suspected Lynch syndrome patients reveals pathogenicity due to attenuated mismatch repair. *J Med Genet*. 51:245-253.
- Ziada-Bouchaar H et al. (2017). First description of mutational analysis of MLH1, MSH2 and MSH6 in Algerian families with suspected Lynch syndrome. *Fam Cancer*. 16:57-66.

P003 Product history

Version	Modification
D1	Two probes specific for the recurrent 10 Mb inversion on chr. 2p have been added; and one reference probe has been replaced.
C1	Two reference probes have been replaced and two extra reference probes have been added.
B2	The 88 and 96 nt DNA denaturation control fragments have been replaced (QDX2).
B1	One MLH1 probe (exon 19) and four MSH2 probes (exons 1, 1, 7 and 10) have been replaced. In addition, one extra MSH2 exon 1 probe, two extra EPCAM (formerly known as TACSTD1) probes and two extra control fragments at 100 and 105 nt have been included.
A2	One extra MSH2 exon 1 probe has been included.
A1	First release.

Implemented changes in the product description

Version D1-04 – 02 August 2018 (04)

- Product description restructured and adapted to a new template.
- Minor textual changes.
- Warning added to Table 2 for probe specificity relying on a single nucleotide difference between target gene and related gene or pseudogene.
- Selected publications were updated using P003.
- Morocco and Israel were added as countries where product has IVD status.

Version D1-03 – 28 November 2017 (03)

- Intended use adjusted to clarify usage of product.
- Text on Binning DNA and Artificial Duplication DNA adjusted.
- Updated references using P003.
- Changed note on the 265 and 317 nt probes.
- Changed note on the 278 nt probe.
- Updated notes under Tables 1 and 2.
- Minor textual changes throughout the document.

Version D1-02 – 10 May 2017 (03)

- Product description restructured and adapted to a new template.

Version D1-01 – 18 November 2015 (02)

- Product description restructured and adapted to a new template.
- Product description adapted to a new product version (version number changed, changes in Table 1 and Table 2).
- Various minor textual or layout changes.

Version 36 (02) – 20 April 2015

- Product description completely rewritten.

Version 35




- Product description adapted to a new lot (lot number added, new picture included).
- Minor textual changes.

Version 34

Product description completely rewritten.

More information: www.mlpa.com; www.mlpa.eu

	MRC-Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands
E-mail	info@mlpa.com (information & technical questions); order@mlpa.com (orders)
Phone	+31 888 657 200

	EUROPE*  MOROCCO ISRAEL
	ALL OTHER COUNTRIES

*comprising EU (candidate) member states and members of the European Free Trade Association (EFTA).
The product is for RUO in all other European countries.