

Product Description SALSA® MLPA® Probemix P088-C2 Oligodendroglioma 1p-19q

To be used with the MLPA General Protocol.

Version C2. As compared to version C1, one reference probe is replaced, lengths of several probes are adjusted with no change in the sequences detected, one reference probes is replaced and sample DNA used for this probemix is changed from SD021 to SD054. For complete product history see page 9.

Catalogue numbers:

- **P088-025R:** SALSA MLPA Probemix P088 Oligodendroglioma 1p-19q, 25 reactions.
- **P088-050R:** SALSA MLPA Probemix P088 Oligodendroglioma 1p-19q, 50 reactions.
- **P088-100R:** SALSA MLPA Probemix P088 Oligodendroglioma 1p-19q, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit, available for various number of reactions. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman capillary sequencers, respectively (see www.mlpa.com).

Certificate of Analysis: Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mlpa.com.

Precautions and warnings: For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mlpa.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

General information: The SALSA MLPA Probemix P088 Oligodendroglioma 1p-19q is a **research use only (RUO)** assay for the detection of deletions and duplications in the 1p and 19q chromosomal regions which are associated with oligodendroglioma. This probemix can also be used to detect the presence of p.R132C and p.R132H point mutations in *IDH1* gene and p.R172K and p.R172M point mutations in *IDH2* gene.

Oligodendrogliomas are central nervous system neoplasms derived from a subset of glial cells, known as oligodendrocytes. The diagnosis of oligodendroglioma is based on the presence of both an *IDH* gene family mutation in combination with co-deletion of 1p and 19q arms (Louis et al. 2016).

This SALSA MLPA Probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.

Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene>

For NM_ mRNA reference sequences: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide>

Locus Reference Genomic (LRG) database: <http://www.lrg-sequence.org/>

Probemix content: The SALSA MLPA Probemix P088-C2 Oligodendroglioma 1p-19q contains 58 MLPA probes with amplification products between 122 and 504 nucleotides (nt). This includes 19 probes for the 1p arm plus three flanking probes for 1q arm, and 11 probes for the 19q arm plus two flanking probes for 19p arm. Furthermore, it also contains four probes specific for the p.R132H and p.R132C mutation of *IDH1* and for the p.R172K and p.R172M mutation of *IDH2* which will only generate a signal when the mutation is present. In addition, 14 reference probes are included and target relatively copy number stable regions in central nervous system tumours, especially in oligodendrogliomas. Complete probe sequences and the identity of the genes detected by the reference probes is available in Table 2b.

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA quantity fragments (Q-fragments), two DNA denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see Table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mlpa.com.

| Length (nt) | Name |
|-------------|--|
| 64-70-76-82 | Q-fragments (Only visible with <100 ng sample DNA) |
| 88-96 | D-fragments (Low signal of 88 or 96 nt fragment indicates incomplete denaturation) |
| 92 | Benchmark fragment |
| 100 | X-fragment (X chromosome specific) |
| 105 | Y-fragment (Y chromosome specific) |

No DNA controls result in only five major peaks shorter than 121 nt: four Q-fragments at 64, 70, 76 and 82 nt, and one 19 nt peak corresponding to the unused portion of the fluorescent PCR primer. Non-specific peaks longer than 121 nt AND with a height >25% of the median of the four Q-fragments should not be observed. Note: peaks below this 25% threshold are not expected to affect MLPA reactions when sufficient amount of sample DNA (50-250 ng) is used.

MLPA technique: The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mlpa.com). More information on the use of MLPA in tumour applications can be found in Hömig-Hölzel and Savola 2012.

Required specimens: Extracted DNA, which includes DNA derived from paraffin-embedded tissues, free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol. More information on the use of FFPE tissue derived DNA for MLPA can be found in Atanesyan et al. 2017.

Reference samples: All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from healthy individuals. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol.

Positive control DNA samples: MRC-Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (<https://catalog.coriell.org>) and DSMZ (<https://www.dsmz.de/home.html>) have a diverse collection of biological resources which may be used as a positive control DNA sample in your MLPA experiments. For example, samples from the Coriell Institute have been tested at MRC-Holland with the P088-C2 probemix; a heterozygous deletion of 1p36.32 (*GNB1*, *TNFRSF14* and *TP73*) and a heterozygous duplication of *CDKN2A* and *CDKN2B* can be detected in NA22976, a heterozygous deletion of 1p36.23 (*PARK7*) can be found in NA50276 and a heterozygous deletion of 1q31.3-1q32.1 (*CRB1* and *TNNT2*) is present in NA00214. The quality of cell lines can change, therefore samples should either be acquired from quality assured biorepositories with minimal cell passage or should be validated before use.

SALSA Binning DNA SD054: The SD054 Binning DNA provided with this probemix can be used as Binning DNA sample for binning of four mutation-specific probes (203 nt probe 19529-L16492 *IDH1* p.R132H=c.395G>A mutation, 227 nt probe 14787-L23353 *IDH1* p.R132C=c.394C>T mutation, 238 nt probe 20963-L29002 *IDH2* p.R172K=c.515G>A mutation and 244 nt probe 20963-L29001 *IDH2* p.R172M=c.515G>T mutation). SD054 Binning DNA is a mixture of genomic DNA from healthy individuals and synthetic DNA that contains the target sequence detected by the above mentioned probes. Inclusion of one reaction with 5 µl SD054 Binning DNA in initial MLPA experiments is essential as it can be used to aid in data binning of the peak pattern using Coffalyser.Net software. Furthermore, Binning DNA should be included in the experiment whenever changes have been applied to the set-up of the capillary electrophoresis device (e.g. when capillaries have been renewed). Binning DNA should never be used as a reference sample in the MLPA data analysis, neither should it be used in quantification of mutation signals, as for this purpose true mutation positive patient samples or cell lines should be used. It is strongly advised to use DNA sample and reference DNA samples extracted with the same method and derived from the same source of tissue. For further details, please consult the SD054 Binning DNA product description provided.

This product is for research use only (RUO).

Data analysis: Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mlpa.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results: The standard deviation of all probes in the reference samples should be ≤ 0.10 . When this criterion is fulfilled, the following cut-off values for the dosage quotient (DQ) of the probes can be used to interpret MLPA results for autosomal or pseudo-autosomal chromosomes:

| Copy Number status | Dosage quotient |
|---|--------------------|
| Normal | $0.80 < DQ < 1.20$ |
| Homozygous deletion | $DQ = 0$ |
| Heterozygous deletion | $0.40 < DQ < 0.65$ |
| Heterozygous duplication | $1.30 < DQ < 1.65$ |
| Heterozygous triplication/ Homozygous duplication | $1.75 < DQ < 2.15$ |
| Ambiguous copy number | All other values |

Please note that these above mentioned dosage quotients are affected both by percentage of tumour cells and by possible subclonality.

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in subclonal cases.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Incomplete DNA denaturation (e.g. due to salt contamination) can lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes are unlikely to have any relation to the condition tested for.

P088 specific note:

- In samples from tumour tissues, reference probes are more prone to have deviating copy number results than in blood derived germline samples. When regions targeted by reference probes are affected by copy number alterations in some cases it can help to turn the slope correction off in Coffalyser.Net analysis to get the correct interpretation of the target region.

Limitations of the procedure:

- In most populations, the most genetic alterations in 1p/19q chromosomal regions are small (point) mutations, most of which will not be detected by using SALSA MLPA Probemix P088 Oligodendroglioma 1p-19q.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.

- Sequence changes (e.g. SNPs, point mutations, small indels) in the target sequence detected by a probe can cause false positive results. Mutations/SNPs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.
- MLPA analysis on tumour samples provides information on the *average* situation in the cells from which the DNA sample was purified. Gains or losses of genomic regions or genes may not be detected if the percentage of tumour cells is low. In addition, subclonality of the aberration affects the final ratio of the corresponding probe. Furthermore, there is always a possibility that one or more reference probes *do* show a copy number alteration in a patient sample, especially in solid tumours with chaotic karyotypes.

Confirmation of results: Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

COSMIC mutation database: <http://cancer.sanger.ac.uk/cosmic>. We strongly encourage users to deposit positive results in the Catalogue Of Somatic Mutations In Cancer (COSMIC). Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report false positive results due to SNPs and unusual results to MRC-Holland: info@mlpa.com.

Table 1. SALSA MLPA Probemix P088-C2 Oligodendroglioma 1p-19q

| Length (nt) | SALSA MLPA probe | Chromosomal position (hg18) | | | | Location (hg18) in kb |
|-------------|--|-----------------------------|----------------|----------|--------------------------|-----------------------|
| | | Reference | Chr. 1 | Chr. 19 | IDH mutation CDKN2A/B | |
| 64-105 | Control fragments – see Table in probemix content section for more information | | | | | |
| 122 | Reference probe 02844-L02274 | 18q11 | | | | 18-019.394 |
| 131 | NOTCH2 probe 05745-L05183 | | 1p12 | | | 01-120.331 |
| 137 | Reference probe 09957-L20646 | 17p13 | | | | 17-007.355 |
| 142 | CDKN2B probe 11867-L23298 | | | | exon 1 | 09-021.999 |
| 148 ↵ | SMARCA4 probe 02488-L22890 | | | 19p13.2 | | 19-011.031 |
| 153 | CDKN2A probe 16881-L23102 | | | | exon 3 | 09-021.961 |
| 157 | UPK1A probe 18116-L23103 | | | 19q13.12 | | 19-040.856 |
| 163 | PTAFR probe 18115-L23104 | | 1p35.3 | | | 01-028.350 |
| 168 « | CCNE1 probe 02881-L23105 | | | 19q12 | | 19-035.005 |
| 173 | Reference probe 15449-L23605 | 12q13 | | | | 12-046.676 |
| 178 | GNB1 probe 02890-L20648 | | 1p36.33 | | | 01-001.747 |
| 184 | PDCD5 probe 02882-L02349 | | | 19q13.11 | | 19-037.764 |
| 190 | CDKN2A probe 16880-L20211 | | | | exon 1 | 09-021.984 |
| 196 ↵ | TNNT2 probe 06557-L20938 | | 1q32.1 | | | 01-199.604 |
| 203 § ¥ | IDH1 probe 19529-L16492 | | | | p.R132H | 02-208.821 |
| 208 | Reference probe 16261-L18553 | 20q11 | | | | 20-034.979 |
| 214 ↵ | LDLR probe 02314-L20213 | | | 19p13.2 | | 19-011.077 |
| 220 | PPP1R15A probe 02887-L02354 | | | 19q13.33 | | 19-054.070 |
| 227 ¥ § | IDH1 probe 14787-L23353 | | | | p.R132C | 02-208.821 |
| 232 | CDKN2B probe 10337-L23606 | | | | exon 2 | 09-021.996 |
| 238 § ¥ | IDH2 probe 20963-L29002 | | | | p.R172K | 15-088.433 |
| 244 § ¥ | IDH2 probe 20963-L29001 | | | | p.R172M | 15-088.433 |
| 253 ¥ | CDKN2A probe 16060-L29163 | | | | exon 2 | 09-021.965 |
| 259 | GTF2B probe 02871-L19715 | | 1p22.2 | | | 01-089.126 |
| 265 Ж « | UPK1A probe 18117-SP0616-L23106 | | | 19q13.12 | | 19-040.861 |
| 270 | Reference probe 16659-L19210 | 2p16 | | | | 02-051.107 |
| 277 * | Reference probe 17450-L29159 | 16p13 | | | | 16-009.761 |
| 283 « | CDKN2C probe 18565-L24220 | | 1p33 | | | 01-051.208 |
| 288 | FAF1 probe 02877-L24219 | | 1p33 | | | 01-051.026 |
| 293 « | WNT4 probe 06055-L24329 | | 1p36.12 | | | 01-022.329 |
| 300 | BAX probe 00348-L00174 | | | 19q13.33 | | 19-054.151 |
| 306 ↵ | LMNA probe 16877-L19710 | | 1q22 | | | 01-154.372 |
| 313 ¥ « | CHMP2A probe 18119-L29136 | | | 19q13.43 | | 19-063.757 |
| 319 ¥ | Reference probe 04833-L22803 | 5p13 | | | | 05-037.032 |
| 326 « | CHMP2A probe 18118-L23300 | | | 19q13.43 | | 19-063.755 |
| 332 « | TP73 probe 01682-L24330 | | 1p36.32 | | | 01-003.558 |
| 340 « | MIR101-1 probe 13654-L24420 | | 1p31.3 | | | 01-065.297 |
| 347 | TNFRSF14 probe 04693-L24421 | | 1p36.32 | | | 01-002.480 |
| 355 | Reference probe 06426-L05952 | 6p22 | | | | 06-024.386 |
| 362 | TGFB1 probe 02889-L23352 | | | 19q13.2 | | 19-046.542 |
| 370 | FUBP1 probe 18571-L24211 | | 1p31.1 | | | 01-078.203 |
| 377 | DPYD probe 02870-L23108 | | 1p21.3 | | | 01-098.159 |
| 385 | Reference probe 08311-L23302 | 11q22 | | | | 11-098.932 |
| 391 | Reference probe 10464-L23212 | 2p11 | | | | 02-085.640 |
| 400 | PRDX1 probe 18410-L23657 | | 1p34.1 | | | 01-045.749 |
| 407 « | CDKN2C probe 18566-L24049 | | 1p33 | | | 01-051.212 |
| 413 ¥ | MFN2 probe 20882-L29180 | | 1p36.22 | | | 01-011.984 |
| 420 | ZNF296 probe 03221-L24213 | | | 19q13.32 | | 19-050.271 |
| 427 ↵ | CRB1 probe 06961-L24214 | | 1q31.3 | | | 01-195.593 |
| 436 | Reference probe 10634-L11182 | 8q12 | | | | 08-061.856 |
| 445 ¥ | PLPP3 probe 18120-L24277 | | 1p32.2 | | | 01-056.775 |
| 454 ¥ | PARK7 probe 02189-L29162 | | 1p36.23 | | | 01-007.968 |
| 463 « | CIC probe 18575-L24215 | | | 19q13.2 | | 19-047.487 |
| 475 | NRAS probe 01032-L20220 | | 1p13.2 | | | 01-115.053 |
| 481 | Reference probe 09772-L10187 | 15q21 | | | | 15-042.706 |
| 489 | Reference probe 17939-L15290 | 3q25 | | | | 03-157.716 |

| Length (nt) | SALSA MLPA probe | Chromosomal position (hg18) | | | Location (hg18) in kb |
|-------------|---------------------------------|-----------------------------|--------|---------|-----------------------|
| | | Reference | Chr. 1 | Chr. 19 | |
| 497 | PRDX1 probe 18413-L24632 | | 1p34.1 | | 01-045.760 |
| 504 | Reference probe 13438-L24633 | 5q31 | | | 05-131.756 |

* New in version C2 (from lot C2-0416 onwards).

‡ Changed in version C2 (from lot C2-0416 onwards). Small change in length, no change in sequence detected.

§ Mutation-specific probe. This probe will only generate a signal when the indicated mutation is present. It has been tested on artificial DNA **but not on positive human samples!**

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

⌘ This probe consists of three parts and has two ligation sites. A low signal of this probe can be due to depurination of the sample DNA, e.g. due to insufficient buffering capacity or a prolonged denaturation time.

→ Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

Note: Please notify us of any mistakes: info@mlpa.com.

Table 2a. P088-C2 target probes arranged according to chromosomal location

| Length (nt) | SALSA MLPA probe | Gene /Exon | Location/ Ligation site | Partial sequence (24 nt adjacent to ligation site) | Distance to next probe | Location (hg18) in kb |
|---|------------------|--|-------------------------|--|------------------------|-----------------------|
| Chromosome 1 | | | | | | |
| Loss of 1p arm is together with loss of 19q arm a diagnostic molecular marker in oligodendrogliomas and it can be used in predicting the response to chemotherapy and prognosis (Smith et al. 2000; Cairncross et al. 1998). | | | | | | |
| 1p arm | | | | | | |
| 178 | 02890-L20648 | GNB1 | 1p36.33 | CTAAGATCGGAA-GATGAGTGAGCT | 733.0 kb | 01-001.747 |
| 347 | 04693-L24421 | TNFRSF14 | 1p36.32 | CAATACCCTCAT-TCACGGGGAGGA | 1.1 Mb | 01-002.480 |
| 332 « | 01682-L24330 | TP73 | 1p36.32 | GAGACCCGGGTG-TCAGGAAAGATG | 4.4 Mb | 01-003.558 |
| 454 | 02189-L29162 | PARK7 | 1p36.23 | AGAGCAGCGAAC-TGCAGCATCAC | 4.0 Mb | 01-007.968 |
| 413 | 20882-L29180 | MFN2 | 1p36.22 | CGCAGAAGGCTT-TCAAGTGAGGAT | 10.3 Mb | 01-011.984 |
| 293 « | 06055-L24329 | WNT4 | 1p36.12 | GCGAGAACTCA-AGGGCCTGATCC | 6.0 Mb | 01-022.329 |
| 163 | 18115-L23104 | PTAFR | 1p35.3 | TGCCCGCTGTA-CCCTTGCAAGAA | 17.4 Mb | 01-028.350 |
| 400 | 18410-L23657 | PRDX1 | 1p34.1 | TACAACCAGTA-GCCTGCCACAA | 10.9 kb | 01-045.749 |
| 497 | 18413-L24632 | PRDX1 | 1p34.1 | ACCTCAGCCATC-CGCAACAGGGTG | 5.3 Mb | 01-045.760 |
| 288 | 02877-L24219 | FAF1 | 1p33 | GGACTGCATTT-AATCCAGCAAGT | 181.6 kb | 01-051.026 |
| 283 « | 18565-L24220 | CDKN2C | 1p33 | CCGGAGGGTTAA-AAGATGATCGCC | 4.3 kb | 01-051.208 |
| 407 « | 18566-L24049 | CDKN2C | 1p33 | TGCTGGAGTTTC-AAGCTGATGTTA | 5.6 Mb | 01-051.212 |
| 445 | 18120-L24277 | PLPP3 | 1p32.2 | AGCACCATCAAG-CCTTACCACCGA | 8.5 Mb | 01-056.775 |
| 340 « | 13654-L24420 | MIR101-1 | 1p31.3 | GGATGGCAGCCA-TCTTACCTTCCA | 12.9 Mb | 01-065.297 |
| 370 | 18571-L24211 | FUBP1 | 1p31.1 | CCATCATGGCGA-TGGACCGGAAA | 10.9 Mb | 01-078.203 |
| 259 | 02871-L19715 | GTF2B | 1p22.2 | CAGATCGGATTT-TAGTGGAGGACT | 9.0 Mb | 01-089.126 |
| 377 | 02870-L23108 | DPYD | 1p21.3 | CTGCTGTCACTT-GGCTCTCTGGCT | 16.9 Mb | 01-098.159 |
| 475 | 01032-L20220 | NRAS | 1p13.2 | TGATGGGACTCA-GGGTTGTATGGG | 5.3 Mb | 01-115.053 |
| 131 | 05745-L05183 | NOTCH2 | 1p12 | GGGGTCAACACT-TACAACCTGCCGC | 34.0 Mb | 01-120.331 |
| 1q arm | | | | | | |
| 306 → | 16877-L19710 | LMNA | 1q22 | ACTGCCTGGCAT-TGTCCAGCTGGA | 41.2 Mb | 01-154.372 |
| 427 → | 06961-L24214 | CRB1 | 1q31.3 | GGAATGTGTGGA-GCTGTCCTCAGA | 4.0 Mb | 01-195.593 |
| 196 → | 06557-L20938 | TNNT2 | 1q32.1 | TTTGCTTCTCT-TCTTCTCATCT | - | 01-199.604 |
| IDH1, at 2q33.3 | | | | | | |
| IDH1 exon numbering according to LRG_610 and indicated ligation sites are in the RefSeq standard NM_005896.3 sequence. The p.R132H (c.395G>A) mutation has been detected by sequencing in 664 samples (92.7%), and the p.R132C (c.394C>T) mutation in 29 samples (4.2%), in a cohort of 1010 WHO grade II and III astrocytomas, oligodendrogliomas and oligoastrocytomas (Hartmann et al. 2009). The IDH1 mutations described are not activating or inactivating but probably result in an altered enzymatic specificity. The probes at 203 nt and 227 nt will only give a signal when respectively the p.R132H or p.R132C mutations are present in the sample. | | | | | | |
| 203 § | 19529-L16492 | IDH1, exon 6; p. R132H (c.395G>A) | 690-691 | CATCATAGGTCA-TCATGCTTATGG | 0.0 kb | 02-208.821 |
| 227 § | 14787-L23353 | IDH1, exon 6; p.R132C (c.394C>T) | 689-688 (reverse) | ATAAGCATGACA-ACCTATGATGAT | - | 02-208.821 |

| Length (nt) | SALSA MLPA probe | Gene /Exon | Location/ Ligation site | Partial sequence (24 nt adjacent to ligation site) | Distance to next probe | Location (hg18) in kb |
|---|---------------------|------------------------------------|-------------------------|--|------------------------|-----------------------|
| CDKN2A/2B, at 9p21.3 | | | | | | |
| CDKN2A exon numbering according to LRG_11 and ligation sites are in the RefSeq standard NM_058195.3, unless stated otherwise. CDKN2B exon numbering and ligation sites as in RefSeq standard NM_078487.2 sequence. Loss of 9p and especially deletions of the 9p21.3 region with CDKN2A/2B are common in oligodendroglioma samples (Bigner et al. 1999) and could be used as a marker of poor prognosis (Cairncross et al. 1998). In contrast, CDKN2A/2B deletions are rare in anaplastic astrocytomas and glioblastomas with mutated IDH1 or IDH2 genes, but are more frequent in these tumours without IDH1/2 mutations (Yan et al. 2009). Additionally, homozygous deletions of CDKN2A/2B have been reported to define a subset of malignant astrocytomas in children (Schiffman et al. 2010). | | | | | | |
| 153 | 16881-L23102 | CDKN2A, exon 3 | 45 nt before exon 3 | TCCTTCCGTC-TGCCGGCCCCCA | 3.7 kb | 09-021.961 |
| 253 | 16060-L29163 | CDKN2A, exon 2 | NM_000077.4; 138-139 | GCCTGGAAGAT-ACCGCGGTCCCT | 19.4 kb | 09-021.965 |
| 190 | 16880-L20211 | CDKN2A, exon 1 | 76-77 | AGTCTGCAGTTA-AGGGGGCAGGAG | 11.4 kb | 09-021.984 |
| 232 | 10337-L23606 | CDKN2B, exon 2 | 1030-1031 | GCCTGTCTGAGA-CTCACAGGAAGG | 3.1 kb | 09-021.996 |
| 142 | 11867-L23298 | CDKN2B, exon 1 | 327-328 | CCAACGGTGGAT-TATCCGGGCCGC | - | 09-021.999 |
| IDH2, at 15q26.1 | | | | | | |
| IDH2 exon numbering according to LRG_611 and ligation sites are in the RefSeq standard NM_002168.3 sequence. The p.R172K (c.515G>A) mutation has been detected by sequencing in 20 samples (64.5%), and the p.R172M (c.515G>T) mutation in 6 samples (19.3%), in a cohort of 1010 WHO grades II and III astrocytomas, oligodendrogliomas and oligoastrocytomas (Hartmann et al. 2009). The same study suggests that IDH2 mutations occur predominantly in oligodendroglial tumours. The IDH2 mutations described are not activating or inactivating but appear to result in an altered enzymatic specificity. The probes at 244 nt and 238 nt will only give a signal when respectively the p.R172M or p.R172K mutations are present in the sample. | | | | | | |
| 238 § | 20963-L29002 | IDH2, exon 5; p.R172K (c. 515 G>A) | 679-680 | TACCATTGGCAA-GCACGCCCATGG | 0.0 kb | 15-088.433 |
| 244 § | 20963-L29001 | IDH2, exon 5; p.R172M (c.515G>T) | 679-680 | TACCATTGGCAT-GCACGCCCATGG | - | 15-088.433 |
| Chromosome 19 | | | | | | |
| Loss of 19q arm is together with loss of 1p arm a diagnostic molecular marker in oligodendrogliomas and it can be used in predicting the response to chemotherapy and prognosis (Smith et al. 2000; Cairncross et al. 1998). The genes FUBP1 (at 1p) and CIC (at 19q) frequently occurred mutated and are suggested to be associated with longer median overall survival (Jiao et al. 2012; Bettgowda et al. 2011). | | | | | | |
| 1p arm | | | | | | |
| 148 ↖ | 02488-L22890 | SMARCA4 | 19p13.2 | CGTCTTGAGTC-GGTCTTACCAG | 45.5 kb | 19-011.031 |
| 214 ↖ | 02314-L20213 | LDLR | 19p13.2 | TCTGTGACTCAG-ACCGGGACTGCT | 23.9 Mb | 19-011.077 |
| 19q arm | | | | | | |
| 168 « | 02881-L23105 | CCNE1 | 19q12 | GATGGTTCATT-TGCCATGGTTAT | 2.8 Mb | 19-035.005 |
| 184 | 02882-L02349 | PDCD5 | 19q13.11 | CGAGGAGCTTGA-GGCCTGAGGAG | 3.1 Mb | 19-037.764 |
| 157 | 18116-L23103 | UPK1A | 19q13.12 | GATGGTGTCAA-CCCATCCCTGAT | 4.8 kb | 19-040.856 |
| 265 Ж « | 18117-SP0616-L23106 | UPK1A | 19q13.12 | GTCTCAGGTGTG-27 nt spanning oligo-ATATCCTTAGCC | 5.7 Mb | 19-040.861 |
| 362 | 02889-L23352 | TGFB1 | 19q13.2 | GAGTGGTTATCT-TTTGATGTCACC | 944.1 kb | 19-046.542 |
| 463 « | 18575-L24215 | CIC | 19q13.2 | GAAACATCCTGC-AGACTGTTGTC | 2.8 Mb | 19-047.487 |
| 420 « | 03221-L24213 | ZNF296 | 19q13.32 | TCATGGACCACA-AGAAGCTGGGCT | 3.8 Mb | 19-050.271 |
| 220 | 02887-L02354 | PPP1R15A | 19q13.33 | GATGTGGATAGT-GAGGATAAGGAA | 81.7 kb | 19-054.070 |
| 300 | 00348-L00174 | BAX | 19q13.33 | TCCCCCGAGAG-GTCTTTTCCGA | 9.6 Mb | 19-054.151 |
| 326 « | 18118-L23300 | CHMP2A | 19q13.43 | TGGAGTTTGAGC-GGCAGGCAGAGA | 2.0 kb | 19-063.755 |
| 313 « | 18119-L29136 | CHMP2A | 19q13.43 | GGGCCCTGAACC-GTCCATGCGGG | - | 19-063.757 |

§ Mutation-specific probe. This probe will only generate a signal when the indicated mutation is present. It has been tested on artificial DNA **but not on positive human samples!**

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

Ж This probe consists of three parts and has two ligation sites. A low signal of this probe can be due to depurination of the sample DNA, e.g. due to insufficient buffering capacity or a prolonged denaturation time.

↖ Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

Note: Exon numbering may differ from literature. When an LRG sequence was available, we have adopted the exon numbering accordingly. Otherwise, the NM sequences indicated in Table 2a are used for exon numbering. Data has been retrieved on 11/2018, but might be changed (e.g. by NCBI) after the release of the product description. Please notify us of any inconsistencies: info@mlpa.com.

Table 2b. Reference probes arranged according to chromosomal location

| Length (nt) | SALSA MLPA probe | Gene | Location | Partial sequence (24 nt adjacent to ligation site) | Location (hg18) in kb |
|-------------|------------------|---------|----------|---|-----------------------|
| 270 | 16659-L19210 | NRXN1 | 2p16 | TAATTTCTGTGG-TTCTTGGGGCTT | 02-051.107 |
| 391 | 10464-L23212 | GGCX | 2p11 | CACCATCATGTT-TCTGGGTGAGGG | 02-085.640 |
| 489 | 17939-L15290 | KCNAB1 | 3q25 | CTTTTCCAGAGA-GAGAAAGTGGAG | 03-157.716 |
| 319 | 04833-L22803 | NIPBL | 5p13 | ACGTGTGAAAAT-GAACAAACGCAA | 05-037.032 |
| 504 | 13438-L24633 | SLC22A5 | 5q31 | GACTTGTATTAT-TTGGCTACAGTC | 05-131.756 |
| 355 | 06426-L05952 | DCDC2 | 6p22 | TTTAGGGAAATG-ATCGCCACTCTA | 06-024.386 |
| 436 | 10634-L11182 | CHD7 | 8q12 | GGATCCCAGTAA-AGGTTTTGGTAA | 08-061.856 |
| 385 | 08311-L23302 | CNTN5 | 11q22 | CACCAGAGCTGT-TAAACACATTGA | 11-098.932 |
| 173 | 15449-L23605 | COL2A1 | 12q13 | CTGGTATCCTCA-TTTTACTTTTTA | 12-046.676 |
| 481 | 09772-L10187 | SPG11 | 15q21 | TTTCTTCAGGAT-TGATAGTCATTC | 15-042.706 |
| 277 | 17450-L29159 | GRIN2A | 16p13 | TGCAGGATTATA-ATCTCACAATCT | 16-009.761 |
| 137 | 09957-L20646 | POLR2A | 17p13 | ACAACAAGAAGA-AGATCATCATCA | 17-007.355 |
| 122 | 02844-L02274 | NPC1 | 18q11 | GACGAGTCTGTG-GATGAGGTCA | 18-019.394 |
| 208 | 16261-L18553 | SAMHD1 | 20q11 | AGTAGACAATGA-GTTGCGTATTTG | 20-034.979 |

Related SALSA MLPA probemixes

- P105 Glioma-2: Contains probes for detection of copy number aberrations of *PDGFRA*, *EGFR*, *CDKN2A*, *PTEN*, *CDK4-MIR26A2-MDM2*, *NFKB1A* and *TP53* genes.
- P370 BRAF-IDH1-IDH2: Contains probes to detect genomic duplications leading to the *KIAA1549-BRAF* and *SRGAP3-RAF1* fusion genes, and identify most common *BRAF*, *IDH1/2* point mutations and copy number alterations of *BRAF*, *CDKN2A/2B*, *FGFR1*, *MYB* and *MYBL1* genes.
- ME012 MGMT-IDH1-IDH2: Contains probes for methylation detection of *MGMT* genes and probes for detection of most common *IDH1/2* point mutations.
- ME024 9p21 CDKN2A/2B region: Contains probes both for methylation and copy number detection for the chromosomal region 9p21 (*CDKN2A/2B*, *MTAP*, *PAX5*).

References

- Atanesyan L et al. (2017). Optimal fixation conditions and DNA extraction methods for MLPA analysis on FFPE tissue-derived DNA. *Am J Clin Pathol.* 147:60-8.
- Bettgowda C et al. (2011). Mutations in CIC and FUBP1 contribute to human oligodendroglioma. *Science.* 333:1453-5.
- Bigner SH et al. (1999). Molecular Genetic Aspects of Oligodendrogliomas Including Analysis by Comparative Genomic Hybridization. *Am J Pathol.* 155:375-86.
- Cairncross JG et al. (1998). Specific genetic predictors of chemotherapeutic response and survival in patients with anaplastic oligodendrogliomas. *J Natl Cancer Inst.* 90:1473-9.
- Hartmann C et al. (2009). Type and frequency of *IDH1* and *IDH2* mutations are related to astrocytic and oligodendroglial differentiation and age: a study of 1,010 diffuse gliomas. *Acta Neuropathol.* 118:469-74.
- Hömig-Hölzel C and Savola S (2012). Multiplex ligation-dependent probe amplification (MLPA) in tumor diagnostics and prognostics. *Diagn Mol Pathol.* 21:189-206.
- Jiao Y et al. (2012). Frequent ATRX, CIC, FUBP1 and *IDH1* mutations refine the classification of malignant gliomas. *Oncotarget.* 3:709-22.
- Louis DN et al. (2016). The 2016 World Health Organization Classification of Tumors of the Central Nervous System: a summary. *Acta Neuropathol.* 131:803-20.
- Schiffman JD et al. (2010). Oncogenic BRAF mutation with *CDKN2A* inactivation is characteristic of a subset of pediatric malignant astrocytomas. *Cancer Res.* 70:512-9.
- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Smith JS et al. (2000). Alterations of chromosome arms 1p and 19q as predictors of survival in oligodendrogliomas, astrocytomas, and mixed oligoastrocytomas. *J Clin Oncol.* 18:636-45.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.
- Yan H et al. (2009). *IDH1* and *IDH2* mutations in gliomas. *N Engl J Med.* 360-765-73.

Selected publications using SALSA MLPA Probemix P088 Oligodendroglioma.

- Bieńkowski M et al. (2018). Molecular diagnostic testing of diffuse gliomas in the real-life setting: A practical approach. *Clin Neuropathol.* 37:166-77.
- Blesa D et al. (2009). Novel genomic alterations and mechanisms associated with tumor progression in oligodendroglioma and mixed oligoastrocytoma. *J Neuropathol Exp Neurol.* 68:274-85.
- Franco-Hernández C et al. (2009). Allelic status of 1p and 19q in oligodendrogliomas and glioblastomas: multiplex ligation-dependent probe amplification versus loss of heterozygosity. *Cancer Genet Cytogenet.* 190:93-6.
- Hench J et al. (2018). Satisfying your neuro-oncologist: a fast approach to routine molecular glioma diagnostics. *Neuro-Oncology.* 20:1682-3.
- Jeuken J et al. (2011). The nature and timing of specific copy number changes in the course of molecular progression in diffuse gliomas: further elucidation of their genetic "life story". *Brain Pathol.* 21:308-20.
- Jeuken J et al. (2006). Multiplex ligation-dependent probe amplification: a diagnostic tool for simultaneous identification of different genetic markers in glial tumors. *J Mol Diagn.* 8:433-43.
- Mellai M et al. (2012). MGMT promoter hypermethylation and its associations with genetic alterations in a series of 350 brain tumors. *J Neurooncol.* 107:617-31.
- Motomura K et al. (2011). Benefits of interferon- β and temozolomide combination therapy for newly diagnosed primary glioblastoma with the unmethylated MGMT promoter: A multicenter study. *Cancer.* 117:1721-30.
- Mukasa A et al. (2012). Significance of IDH mutations varies with tumor histology, grade, and genetics in Japanese glioma patients. *Cancer Sci.* 103:587-92.
- Natte R et al. (2005). Multiplex ligation-dependent probe amplification for the detection of 1p and 19q chromosomal loss in oligodendroglial tumors. *Brain Pathol.* 15:192-7.
- Shibahara I et al. (2012). *IDH1/2* gene status defines the prognosis and molecular profiles in patients with grade III gliomas. *Int J Clin Oncol.* 17:551-61.
- Weller M et al. (2007). Combined 1p/19q loss in oligodendroglial tumors: predictive or prognostic biomarker? *Clin Cancer Res.* 13:6933-7.

| P088 Product history | |
|-----------------------------|---|
| <i>Version</i> | <i>Modification</i> |
| <i>C2</i> | <i>Lengths of several probes are adjusted with no change in the sequences detected, one reference probes is replaced and sample DNA used for this probemix is changed from SD021 to SD054.</i> |
| <i>C1</i> | <i>Several target probes and all reference probes are replaced. Also, four probes for point mutations in IDH1 and IDH2, probes for CDKN2A/2B and new control fragments (QDX2) are included.</i> |
| <i>B2</i> | <i>The 88 and 96 nt control fragments have been replaced (QDX2).</i> |
| <i>B1</i> | <i>Several probes have been replaced, incl. three probes for 1q. Extra control fragments have been included.</i> |
| <i>A1</i> | <i>First release.</i> |

Implemented changes in the product description
Version C2-01 – 04 December 2018 (01P)

- Product description restructured and adapted to a new template.
- For uniformity, the chromosomal positions and bands in this document are now all based on hg18 (NCBI36).
- Warning about SNPs removed in Table 1 and 2 for NOTCH2 (05745-L05183), CDKN2B (10337-L23606) and ZNF296 (03221-L24213) probes.

Version 21 – 15 August 2017 (T08)

- One new reference for P088 added on page 1.
- Small layout changes and corrections of typos in the document.


Version 20 – 2 May 2016 (T08)

- Product description adapted to a new product version (version number changed, lot number added, small changes in Table 1 and Table 2, new picture included).
- A new SD is available for P088-C2-0416. The previous SD021 is changed to SD054.
- New references added on pages 1 and 2.
- Warning added in Table 1, SNP rs531705888 could influence the probe signal at probe 232nt 10337-L23606.
- PLPP3 gene name has been changed (PPAP2B previously) in Table 1 and 2.

Version 19 – 22 May 2015 (T07)

- Information on the p.R172M (G515T) mutation corrected in Table 2.
- Minor restructuring of the document.

More information: www.mlpa.com; www.mlpa.eu

| | |
|---|---|
|  | MRC-Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands |
| E-mail | info@mlpa.com (information & technical questions); order@mlpa.com (orders) |
| Phone | +31 888 657 200 |